

Product Introduction

The DV-300 Aspirator is a high quality Portable Aspirator intended for use in accordance with the instructions of a healthcare professional or physician. The device is designed with a maintenance-free pump and adjustable vacuum settings. The DV-300 is a high vacuum/ low flow unit.

Intended Use / Indications for Use

The device is to be used to remove fluids from the treated tracheotomy patient airway.

Device Description

The DV-300 Aspirator is intended to be used to remove fluids from the treated tracheotomy patient airway. It creates a negative pressure (vacuum) that draws the fluids through disposable tubing that is connected to a collection container. The fluids are trapped in the collection container for the proper disposal. The device is only for use on the order of a physician. If practiced outside of the hospital environment, the care giver must be trained by a certified healthcare professional and the training must be recorded.

The DV-300 Aspirator is a portable AC-powered suction pump. It consists of an on/off switch, a pump unit, a power cord, a collection container, relief valve, pressure gauge, pressure adjustment knob, inline filter, and suction tubing.

Product Specification

Power requirement: AC 110-120V 60Hz

Max working air flow: ~18 l/min

Vacuum: 0 - 620 mmHg

Dimensions: L34 x W16.5 x H19 cm

Weight: 3.5 kg / 7.7 lbs

Collection container capacity: 800 ml

Power cord length: 180 cm

Operating temp: 0 - 40 °C

Operating relative humidity: 0 - 90%

Storage and transport temp: -20 - 50 °C


Storage and transport humidity: 0 - 95%

Operating atmospheric pressure: 10.2 - 15.4 psi (70 - 106 kPa)

Storage atmospheric pressure: 7.3 - 15.4 psi (50 - 106 kPa)

Part List

- | | | |
|--------------------------------|---------------------------------|---------------------------|
| 1. Intermediate tubing | 6. Collection container | 11. Handle |
| 2. Inline filter | 7. Collection container coaster | 12. Power cord and plug * |
| 3. Inlet | 8. Power switch | 13. Power inlet |
| 4. Collection container holder | 9. Vacuum gauge | 14. Vacuum port |
| 5. Suction tubing | 10. Vacuum regulator knob | 15. Patient port |



* Illustration of plug is for reference only. Plug adaptor may vary depending on your region.