



AURACAST RESOURCE CENTER FROM OAKTREE PRODUCTS

Auracast Pairing Guide

App navigation for compatible hearing devices

Use this quick reference to locate Auracast options in hearing device apps. Compatibility and menu names vary by model, firmware, app version, and smartphone. Confirm details with the manufacturer or hearing care provider.

Brand and model	Joining method	App navigation pathway
 life is on EON / Infinio	QR Code Only via App	myPhonak app > Program Menu > "Join an Auracast broadcast"
 Nucleus 8 / Baha 7	QR Code Only via App	Nucleus Smart app > Audio Sharing > Scan QR Code
 Nexia / Sensia / Vivia	Live App Text List	ReSound Smart 3D app > More > "Find a broadcast"
 Intent / Zeal / Reveal	Live App Text List	Oticon Companion app > More > "Auracast broadcast"
 Edge AI / Omega AI	Live App Text List	My Starkey app > Control > Audio Streams > Auracast
 Pure Charge&Go IX	Live App Text List	Signia app > Settings > Connectivity > Auracast Stream
Classic / Older Hearing Aids	Listen Tech RX1 + Neck Loop	No Smartphone Required (Switch hearing aid to "T" position)

Listening tips

Volume: Use the broadcast volume controls available in the hearing device app or on the hearing devices. Phone volume-button behavior can vary by setup.

Android: On supported phones and hearing devices, system Auracast controls or a broadcast QR code may offer another joining route. Android version alone does not guarantee compatibility.

Logos and product names are trademarks of their respective owners.