



How to Join an Auracast Broadcast

Choose the connection method that works with your listening device.

01 Use the Auri receiver

The AURI-RX1 receiver connects to an Auracast broadcast. The LA-470 neck loop can be used with telecoil-enabled hearing devices.

1. If scanning is available, press the right button below the magnifying glass to scan.
2. Use the volume buttons or left button to find your broadcast.
3. Press the right button to select the broadcast.



02 Scan a venue QR code

Look for printed QR codes near your seating, gate, or service area.

1. Unlock your phone and scan the venue's QR code with your phone camera.
2. Follow the prompt on your phone to join the broadcast.



Scan to listen to the Hockey channel

03 Use your phone's settings

ANDROID 15+

1. Make sure compatible earbuds or hearing aids are on and connected.
2. Go to Settings → Bluetooth.
3. Tap the gear next to your earbuds' name.
4. Look for "Find a broadcast," "Listen to Auracast," or "Audio streams."
5. Select the stream name to start listening.



04 Use your hearing device's app

1. Open the official app you normally use for your hearing aids or earbuds (e.g., ReSound, Oticon, or Philips).
2. Look for Auracast, Audio Sharing, or Broadcast Search.
3. Tap it to search for the audio stream and listen.



Steps and menu names may vary by phone, app, and listening device.